Restore Foam Dressing without Border

with TRIACT ADVANCED Technology

DESCRIPTION

Restore Foam Dressing without Border is a non-adhesive absorbent, micro-adherent, hydrocellular dressing featuring TRIACT ADVANCED Technology. Restore Foam Dressing without Border consists of:

- A pre-formed grid containing a micro-adherent lipido-colloid matrix comprised of hydrocolloid particles (carboxymethylcellulose) and cohesion polymers that remain in contact with the wound.
- The lipido-colloid matrix is coated on a non-sensitizing, super-absorbent polyurethane foam pad.
- The foam pad has a non-woven protective backing made of elastic polyurethane.

Restore Foam Dressing without Border is supplied in sterile, individual packaging for single use.

INDICATIONS FOR USE

Restore Foam Dressing without Border is indicated in the treatment of moderate to heavily exuding wounds, including chronic wounds (pressure sores/ulcers, diabetic ulcers, venous stasis ulcers, arterial ulcers, leg ulcers of mixed etiology, malignant wounds, etc.) and acute wounds (traumatic wounds, surgical wounds, donor sites, partial-thickness burns, skin abrasions, lacerations, minor cuts, scalds, etc.). Due to its non-adhesive nature, Restore Foam Dressing without Border is recommended in the treatment of wounds in which the surrounding skin is fragile, friable or damaged. It is suitable for use under compression bandages / wraps due to the ability of the dressing to retain exudate.

MECHANISM OF ACTION

In the presence of exudate, the hydrocolloid particles in the TRIACT ADVANCED Technology polymer matrix form a lipido-colloid gel at the interface between the wound and the dressing. This results in:

- Virtually pain-free dressing changes for the patient
- Maintenance of a moist environment that promotes healing
- Gel formation in contact with the wound
- Minimizes trauma, such as pain and bleeding, which may be experienced during dressing removal

The backing is soft, flexible and very conformable to allow the dressing to be easily shaped to the anatomical contours of the wound.

The super-absorbent foam pad allows containment of exudate and protects the periwound skin from maceration.

INSTRUCTIONS FOR USE

1. Wound preparation:
 - Cleanse the wound in accordance with established procedures.
 - Carefully dry the periwound tissue.
2. Dressing application:
 - Select the appropriate dressing size to ensure the dressing will cover the entire wound.
 - The Restore Foam Dressing without Border can be cut with the release liner in place in order to optimize placement over wounds on heels, elbows, etc.
 - Remove the protective tabs from the dressing.
 - Apply the micro-adherent side of the dressing to the wound.
 - Secure the dressing using a bandage or a tape.
 - The dressing can be applied under compression.
3. Dressing changes:
 - Restore Foam Dressing without Border can be left in place up to 7 days, depending on the exudate and the clinical condition of the wound.
 - Duration of treatment is determined by the physician and depends on wound type and conditions.

PRECAUTIONS FOR USE

- Do not re-use the dressing.
- The micro-adherent mass of the Restore Foam Dressing without Border pad sticks to latex surgical gloves. It is therefore recommended to use latex-free gloves.
- Do not re-sterilize the dressing.
- Store Restore Foam Dressing without Border flat, away from moisture and heat (temperature below 25°C/77°F). Protect from light.

CONTRAINDICATIONS

Restore Foam Dressing without Border should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or one of its components.

HOW SUPPLIED

REF.: 509380: 2.5" x 2.5" (6 cm x 6 cm)
509381: 4" x 4" (10 cm x 10 cm)
509382: 6" x 6" (15 cm x 15 cm)
509383: 8" x 8" (20 cm x 20 cm)
509384: Heel 4.7" x 7.5" (12 cm x 19 cm)

Each dressing is individually packed in a sterile pouch. Sterilization is by radiation. Sterility is guaranteed unless a package is damaged or opened. Single use only.

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Latest revision : 01/2010

hollisterwoundcare

Manufactured for:
Hollister Wound Care, LLC
Libertyville, Illinois 60048 USA

Distributed in Canada by:
Distribué au Canada par:
Hollister Limited
Aurora, Ontario L4G 1G3

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©2010 Hollister Incorporated
Hollisterwoundcare and wave logo are trademarks of Hollister Incorporated.
Restore, TRIACT ADVANCED and graphic are trademarks of Hollister Wound Care, LLC.